



TLR9 Antibody, Biotin conjugated

Product Code	CSB-PA023608LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR96
Immunogen	Recombinant Human Toll-like receptor 9 protein (425-563AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 9 (CD antigen CD289), TLR9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TLR9