



WNT1 Antibody, FITC conjugated

Product Code	CSB-PA026128LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04628
Immunogen	Recombinant Human Proto-oncogene Wnt-1 protein (142-288AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proto-oncogene Wnt-1 (Proto-oncogene Int-1 homolog), WNT1, INT1
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	WNT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.