



DNAJB6 Antibody, Biotin conjugated

Product Code	CSB-PA03239D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75190
Immunogen	Recombinant Human DnaJ homolog subfamily B member 6 protein (208-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily B member 6 (HHDJ1) (Heat shock protein J2) (HSJ-2) (MRJ) (MSJ-1), DNAJB6, HSJ2 MRJ MSJ1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DNAJB6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.