



LEPR Antibody, Biotin conjugated

Product Code	CSB-PA06809D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48357
Immunogen	Recombinant Human Leptin receptor protein (874-1075AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leptin receptor (LEP-R) (HuB219) (OB receptor) (OB-R) (CD antigen CD295), LEPR, DB OBR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LEPR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.