



ICAM1 Antibody, HRP conjugated

Product Code	CSB-PA07149B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05362
Immunogen	Recombinant Human Intercellular adhesion molecule 1 protein (320-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intercellular adhesion molecule 1 (ICAM-1) (Major group rhinovirus receptor) (CD antigen CD54), ICAM1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ICAM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.