




MMP8 Antibody

Product Code	CSB-PA07469A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22894
Immunogen	Recombinant Human Neutrophil collagenase protein (22-283AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:10000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil collagenase (EC 3.4.24.34) (Matrix metalloproteinase-8) (MMP-8) (PMNL collagenase) (PMNL-CL), MMP8, CLG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP8

Image



Western Blot

Positive WB detected in: Mouse liver tissue,
Mouse lung tissue, Mouse kidney tissue


All lanes: MMP8 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 54 kDa

Observed band size: 54 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA07469A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.