



MMP8 Antibody, HRP conjugated

Product Code	CSB-PA07469B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22894
Immunogen	Recombinant Human Neutrophil collagenase protein (22-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil collagenase (EC 3.4.24.34) (Matrix metalloproteinase-8) (MMP-8) (PMNL collagenase) (PMNL-CL), MMP8, CLG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.