



KL Antibody, Biotin conjugated

Product Code	CSB-PA12729D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UEF7
Immunogen	Recombinant Human Klotho protein (34-981AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Klotho (EC 3.2.1.31) [Cleaved into: Klotho peptide], KL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.