




NOS1 Antibody

Product Code	CSB-PA14099A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29475
Immunogen	Recombinant Human Nitric oxide synthase, brain protein (655-875AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:10000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitric oxide synthase, brain (EC 1.14.13.39) (Constitutive NOS) (NC-NOS) (NOS type I) (Neuronal NOS) (N-NOS) (nNOS) (Peptidyl-cysteine S-nitrosylase NOS1) (bNOS), NOS1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOS1

Image



Western Blot

Positive WB detected in: Mouse spleen tissue,

Mouse lung tissue, Mouse kidney tissue


All lanes: NOS1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 161, 149, 126, 44, 165 kDa

Observed band size: 161, 126, 92 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA14099A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.