



FOXO1 Antibody, Biotin conjugated

Product Code	CSB-PA16019D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12778
Immunogen	Recombinant Human Forkhead box protein O1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein O1 (Forkhead box protein O1A) (Forkhead in rhabdomyosarcoma), FOXO1, FKHR FOXO1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FOXO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.