





AQP6 Antibody

Product Code	CSB-PA16139A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13520
Immunogen	Recombinant Human Aquaporin-6 protein (1-30AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-6 (AQP-6) (Aquaporin-2-like) (Kidney-specific aquaporin) (hKID), AQP6, AQP2L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP6

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA16139A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA16139A0Rb at dilution of 1:100