




MYBPC3 Antibody

Product Code	CSB-PA613524HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14896
Immunogen	Recombinant Human Myosin-binding protein C, cardiac-type protein (1-328AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-binding protein C, cardiac-type (Cardiac MyBP-C) (C-protein, cardiac muscle isoform), MYBPC3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MYBPC3

Image



Western Blot

Positive WB detected in: Mouse heart tissue

All lanes: MYBPC-3 antibody at 2.5μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 141 kDa

Observed band size: 141 kDa