



TRIM25 Antibody, Biotin conjugated

Product Code	CSB-PA617909LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14258
Immunogen	Recombinant Human E3 ubiquitin/ISG15 ligase TRIM25 protein (81-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin/ISG15 ligase TRIM25 (EC 6.3.2.n3) (Estrogen-responsive finger protein) (RING finger protein 147) (RING-type E3 ubiquitin transferase) (EC 2.3.2.27) (RING-type E3 ubiquitin transferase TRIM25) (Tripartite motif-containing protein 25) (Ubiquitin/ISG15-conjugating enzyme TRIM25) (Zinc finger protein 147), TRIM25, EFP RNF147 ZNF147
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.