



DOC2B Antibody, Biotin conjugated

Product Code	CSB-PA618877LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14184
Immunogen	Recombinant Human Double C2-like domain-containing protein beta protein (7-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Double C2-like domain-containing protein beta (Doc2-beta), DOC2B, DOC2BL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DOC2B