



DPYD Antibody, HRP conjugated

Product Code	CSB-PA622516LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12882
Immunogen	Recombinant Human Dihydropyrimidine dehydrogenase [NADP(+)] protein (5-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropyrimidine dehydrogenase [NADP(+)] (DHPDHase) (DPD) (EC 1.3.1.2) (Dihydrothymine dehydrogenase) (Dihydrouracil dehydrogenase), DPYD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPYD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.