



MTAP Antibody, FITC conjugated

Product Code	CSB-PA622639LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13126
Immunogen	Recombinant Human S-methyl-5'-thioadenosine phosphorylase protein (1-283AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S-methyl-5'-thioadenosine phosphorylase (EC 2.4.2.28) (5'-methylthioadenosine phosphorylase) (MTA phosphorylase) (MTAP) (MTAPase), MTAP, MSAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTAP