



Plbd2 Antibody, HRP conjugated

Product Code	CSB-PA683927LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4QQW8
Immunogen	Recombinant Rat Putative phospholipase B-like 2 protein (36-585AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative phospholipase B-like 2 (EC 3.1.1.-) (LAMA-like protein 2) (Lamina ancestor homolog 2) (Phospholipase B domain-containing protein 2), Plbd2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Plbd2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.