



PLEKHO1 Antibody, Biotin conjugated

Product Code	CSB-PA706633LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53GL0
Immunogen	Recombinant Human Pleckstrin homology domain-containing family O member 1 protein (143-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family O member 1 (PH domain-containing family O member 1) (C-Jun-binding protein) (JBP) (Casein kinase 2-interacting protein 1) (CK2-interacting protein 1) (CKIP-1) (Osteoclast maturation-associated gene 120 protein), PLEKHO1, CKIP1 OC120
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLEKHO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.