



TDRD12 Antibody, FITC conjugated

Product Code	CSB-PA707440LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q587J7
Immunogen	Recombinant Human Putative ATP-dependent RNA helicase TDRD12 protein (991-1164AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative ATP-dependent RNA helicase TDRD12 (EC 3.6.4.13) (ES cell-associated transcript 8 protein) (Tudor domain-containing protein 12), TDRD12, ECAT8
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TDRD12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.