



IL23R Antibody, Biotin conjugated

Product Code	CSB-PA733162LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VWK5
Immunogen	Recombinant Human Interleukin-23 receptor protein (15-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-23 receptor (IL-23 receptor) (IL-23R), IL23R
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL23R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.