



NOTUM Antibody, HRP conjugated

Product Code	CSB-PA747435LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P988
Immunogen	Recombinant Human Palmitoleoyl-protein carboxylesterase NOTUM protein (100-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoleoyl-protein carboxylesterase NOTUM (EC 3.1.1.98) (hNOTUM), NOTUM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOTUM