




GRAMD4 Antibody

Product Code	CSB-PA764225LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IC98
Immunogen	Recombinant Human GRAM domain-containing protein 4 protein (23-161AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GRAM domain-containing protein 4 (Death-inducing protein), GRAMD4, DIP KIAA0767
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GRAMD4

Image



Western Blot

Positive WB detected in: Mouse brain tissue


All lanes: GRAMD4 antibody at 3.4µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67, 11 kDa

Observed band size: 67 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA764225LA01HU at dilution of 1:100



Western Blot

Positive WB detected in: MCF-7 whole cell lysate

All lanes: GRAMD4 antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67, 11 kDa

Observed band size: 67 kDa