



GRAMD4 Antibody, Biotin conjugated

Product Code	CSB-PA764225LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IC98
Immunogen	Recombinant Human GRAM domain-containing protein 4 protein (23-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GRAM domain-containing protein 4 (Death-inducing protein), GRAMD4, DIP KIAA0767
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GRAMD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.