



FILIP1 Antibody, Biotin conjugated

Product Code	CSB-PA800226LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7B0
Immunogen	Recombinant Human Filamin-A-interacting protein 1 protein (957-1192AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Filamin-A-interacting protein 1 (FILIP), FILIP1, KIAA1275
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FILIP1