



PHLDB2 Antibody, Biotin conjugated

Product Code	CSB-PA803110LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SQ0
Immunogen	Recombinant Human Pleckstrin homology-like domain family B member 2 protein (734-1030AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology-like domain family B member 2 (Protein LL5-beta), PHLDB2, LL5B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PHLDB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.