



SLC37A2 Antibody, Biotin conjugated

Product Code	CSB-PA819478LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TED4
Immunogen	Recombinant Human Glucose-6-phosphate exchanger SLC37A2 protein (46-84AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphate exchanger SLC37A2 (Solute carrier family 37 member 2), SLC37A2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC37A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.