



PHIP Antibody

Product Code	CSB-PA837866LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWQ0
Immunogen	Recombinant Human PH-interacting protein (508-608AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PH-interacting protein (PHIP) (IRS-1 PH domain-binding protein) (WD repeat-containing protein 11), PHIP, WDR11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PHIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.