



TRAPPC9 Antibody, FITC conjugated

Product Code	CSB-PA842741LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q05
Immunogen	Recombinant Human Trafficking protein particle complex subunit 9 protein (780-1030AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking protein particle complex subunit 9 (NIK- and IKBKB-binding protein) (Tularik gene 1 protein), TRAPPC9, KIAA1882 NIBP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRAPPC9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.