



SPPL3 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA848402LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8TCT6 |
| Immunogen | Recombinant Human Signal peptide peptidase-like 3 protein (213-263AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Signal peptide peptidase-like 3 (SPP-like 3) (EC 3.4.23.-) (Intramembrane protease 2) (IMP-2) (Presenilin homologous protein 1) (PSH1) (Presenilin-like protein 4), SPPL3, IMP2 PSL4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | SPPL3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |