



SPPL3 Antibody, Biotin conjugated

Product Code	CSB-PA848402LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCT6
Immunogen	Recombinant Human Signal peptide peptidase-like 3 protein (213-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptide peptidase-like 3 (SPP-like 3) (EC 3.4.23.-) (Intramembrane protease 2) (IMP-2) (Presenilin homologous protein 1) (PSH1) (Presenilin-like protein 4), SPPL3, IMP2 PSL4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPPL3