




# DGAT2 Antibody

<b>Product Code</b>	CSB-PA850407LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96PD7
<b>Immunogen</b>	Recombinant Human Diacylglycerol O-acyltransferase 2 protein (113-388AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	DGAT2

## Image



### Western Blot

Positive WB detected in: HepG2 whole cell lysate, U251 whole cell lysate, A549 whole cell lysate, HeLa whole cell lysate, MCF7 whole cell lysate, K562 whole cell lysate

All lanes: DGAT2 antibody at 1:1000

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa