



DNAAF4 Antibody, Biotin conjugated

Product Code	CSB-PA855504LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXU2
Immunogen	Recombinant Human Dynein assembly factor 4, axonemal protein (109-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein assembly factor 4, axonemal (Dyslexia susceptibility 1 candidate gene 1 protein), DNAAF4, DYX1C1 EKN1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	DNAAF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.