



ZFP42 Antibody, Biotin conjugated

Product Code	CSB-PA856985LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MM3
Immunogen	Recombinant Human Zinc finger protein 42 homolog protein (6-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 42 homolog (Zfp-42) (Reduced expression protein 1) (REX-1) (hREX-1) (Zinc finger protein 754), ZFP42, REX1 ZNF754
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	ZFP42