



ITSN2 Antibody, Biotin conjugated

Product Code	CSB-PA873709LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZM3
Immunogen	Recombinant Human Intersectin-2 protein (940-1209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intersectin-2 (SH3 domain-containing protein 1B) (SH3P18) (SH3P18-like WASP-associated protein), ITSN2, KIAA1256 SH3D1B SWAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITSN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.