



# ITSN2 Antibody, Biotin conjugated

Product Code	CSB-PA873709LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZM3
Immunogen	Recombinant Human Intersectin-2 protein (940-1209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intersectin-2 (SH3 domain-containing protein 1B) (SH3P18) (SH3P18-like WASP-associated protein), ITSN2, KIAA1256 SH3D1B SWAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITSN2