



EHD2 Antibody, Biotin conjugated

Product Code	CSB-PA873710LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN4
Immunogen	Recombinant Human EH domain-containing protein 2 protein (401-543AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EH domain-containing protein 2 (PAST homolog 2), EHD2, PAST2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EHD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.