



MS4A6A Antibody, HRP conjugated

Product Code	CSB-PA875659LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2W1
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 6A protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 6A (CD20 antigen-like 3) (Four-span transmembrane protein 3), MS4A6A, 4SPAN3 CD20L3 MS4A6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MS4A6A