



MS4A6A Antibody, Biotin conjugated

Product Code	CSB-PA875659LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2W1
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 6A protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 6A (CD20 antigen-like 3) (Four-span transmembrane protein 3), MS4A6A, 4SPAN3 CD20L3 MS4A6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MS4A6A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.