



FZD8 Antibody, Biotin conjugated

Product Code	CSB-PA880970LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H461
Immunogen	Recombinant Human Frizzled-8 protein (55-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Frizzled-8 (Fz-8) (hFz8), FZD8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FZD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.