





SLCO1B3 Antibody

Product Code	CSB-PA882065LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPD5
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 1B3 protein (648-695AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 1B3 (Liver-specific organic anion transporter 2) (LST-2) (Organic anion transporter 8) (Organic anion-transporting polypeptide 8) (OATP-8) (Solute carrier family 21 member 8), SLCO1B3, LST2 OATP1B3 OATP8 SLC21A8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLCO1B3


Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA882065LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA882065LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA882065LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)