



SLCO1B3 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA882065LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NPD5 |
| Immunogen | Recombinant Human Solute carrier organic anion transporter family member 1B3 protein (648-695AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Solute carrier organic anion transporter family member 1B3 (Liver-specific organic anion transporter 2) (LST-2) (Organic anion transporter 8) (Organic anion-transporting polypeptide 8) (OATP-8) (Solute carrier family 21 member 8), SLCO1B3, LST2 OATP1B3 OATP8 SLC21A8 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | SLCO1B3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |