



SLCO1B3 Antibody, Biotin conjugated

Product Code	CSB-PA882065LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPD5
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 1B3 protein (648-695AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 1B3 (Liver-specific organic anion transporter 2) (LST-2) (Organic anion transporter 8) (Organic anion-transporting polypeptide 8) (OATP-8) (Solute carrier family 21 member 8), SLCO1B3, LST2 OATP1B3 OATP8 SLC21A8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLCO1B3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.