



USP25 Antibody, HRP conjugated

Product Code	CSB-PA883400LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHP3
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 25 protein (348-606AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 25 (EC 3.4.19.12) (Deubiquitinating enzyme 25) (USP on chromosome 21) (Ubiquitin thioesterase 25) (Ubiquitin-specific-processing protease 25), USP25, USP21
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.