



SLCO3A1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA883432LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UIG8 |
| Immunogen | Recombinant Human Solute carrier organic anion transporter family member 3A1 protein (439-539AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Solute carrier organic anion transporter family member 3A1 (OATP3A1) (Organic anion transporter polypeptide-related protein 3) (OATP-RP3) (OATPRP3) (Organic anion-transporting polypeptide D) (OATP-D) (PGE1 transporter) (Sodium-independent organic anion transporter D) (Solute carrier family 21 member 11), SLCO3A1, OATP3A1 OATPD SLC21A11 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | SLCO3A1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |