



SLCO3A1 Antibody, Biotin conjugated

Product Code	CSB-PA883432LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIG8
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 3A1 protein (439-539AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 3A1 (OATP3A1) (Organic anion transporter polypeptide-related protein 3) (OATP-RP3) (OATPRP3) (Organic anion-transporting polypeptide D) (OATP-D) (PGE1 transporter) (Sodium-independent organic anion transporter D) (Solute carrier family 21 member 11), SLCO3A1, OATP3A1 OATPD SLC21A11
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLCO3A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.