



SLCO3A1 Antibody, Biotin conjugated

Product Code	CSB-PA883432LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIG8
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 3A1 protein (439-539AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 3A1 (OATP3A1) (Organic anion transporter polypeptide-related protein 3) (OATP-RP3) (OATPRP3) (Organic anion-transporting polypeptide D) (OATP-D) (PGE1 transporter) (Sodium-independent organic anion transporter D) (Solute carrier family 21 member 11), SLCO3A1, OATP3A1 OATPD SLC21A11
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLCO3A1