



SIRT7 Antibody, FITC conjugated

Product Code	CSB-PA885703LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRC8
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-7 protein (242-379AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase sirtuin-7 (EC 3.5.1.-) (Regulatory protein SIR2 homolog 7) (SIR2-like protein 7), SIRT7, SIR2L7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SIRT7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.