



# SENP1 Antibody, Biotin conjugated

|                     |                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------|
| Product Code        | CSB-PA885787LD01HU                                                                      |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                           |
| Uniprot No.         | Q9P0U3                                                                                  |
| Immunogen           | Recombinant Human Sentrin-specific protease 1 protein (63-214AA)                        |
| Raised In           | Rabbit                                                                                  |
| Species Reactivity  | Human                                                                                   |
| Tested Applications | ELISA                                                                                   |
| Form                | Liquid                                                                                  |
| Conjugate           | Biotin                                                                                  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4        |
| Purification Method | >95%, Protein G purified                                                                |
| Isotype             | IgG                                                                                     |
| Clonality           | Polyclonal                                                                              |
| Alias               | Sentrin-specific protease 1 (EC 3.4.22.-) (Sentrin/SUMO-specific protease SENP1), SENP1 |
| Immunogen Species   | Homo sapiens (Human)                                                                    |
| Research Area       | Cell Biology                                                                            |
| Target Names        | SENP1                                                                                   |