



FBXO9 Antibody, FITC conjugated

Product Code	CSB-PA891719LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK97
Immunogen	Recombinant Human F-box only protein 9 protein (13-111AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 9 (Cross-immune reaction antigen 1) (Renal carcinoma antigen NY-REN-57), FBXO9, FBX9 VCIA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.