




AJAP1 Antibody

Product Code	CSB-PA891947LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKB5
Immunogen	Recombinant Human Adherens junction-associated protein 1 protein (308-394AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adherens junction-associated protein 1 (Membrane protein shrew-1), AJAP1, MOT8 SHREW1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AJAP1

Image



Western Blot

Positive WB detected in: Mouse brain tissue,
Mouse skeletal muscle tissue


All lanes: AJAP1 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA891947LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.